


TC30

Umgekehrter Kegel • Inverted Cone
Cône renversé

max. 300.000

5



Shank	L mm						
Size	ISO	Order No.	0,7	0,8	0,9	1,0	1,2
Ø ₁₀ mm			006	008	009	010	012
FG	500 314 010 175...	TC30.314...	006	008	009	010	012
							014




TC41

Rund • Round • Rond


- 12 Blades - 010/012/014/016
- 18 Blades - 018
- 20 Blades - 023/027
- 24 Blades - 031
- 30 Blades - 018

FG 300.000, RA 100.000
as from Ø 023, FG 140.000,
RA 50.000

5



Shank	L mm								
Size	ISO	Order No.	0,8	1,0	1,2	1,4	1,6	2,06	2,45
Ø ₁₀ mm		US No.	010	012	014	016	018	023	027
FG	• 500 314 001 071...	TC41.314...	7002	7003	7004		7006	7008	7009
	• 500 314 001 031...	TC41UF.314...	010	012	014	016	018	023	027
RA	• 500 204 001 071...	TC41.204...	010	012	014	016	018	023	027
	• 500 204 001 031...	TC41UF.204...					018		031




TC46

Flamme • Flame • Flamme
• 12 Schneiden • Blades • Lames

FG 300.000, RA 100.000

5



Shank	L mm				
Size	ISO	Order No.	3,5	3,5	3,8
Ø ₁₀ mm		US No.	012	014	018
FG	• 500 314 254 072...	TC46.314...	7103	7104	7106
RA	• 500 204 254 072...	TC46.204...	012	014	018
			012	014	023

TC47

Birne • Pear • Poire
• 12 Schneiden • Blades • Lames

300.000

5




Shank	L mm		1,7	1,8	2,4
Size	ISO	Order No.	010	014	018
Ø ₁₀ mm					
FG	• 500 314 233 072...	TC47.314...	010	014	018

TC47L

Birne lang • Long Pear • Poire long
• 12 Schneiden • Blades • Lames

max. 300.000

5



Shank	L mm				
Size	ISO	Order No.	3,8	4,2	4,2
Ø ₁₀ mm		US No.	010	012	014
FG	• 500 314 234 072...	TC47L.314...	7302	7303	7304
			010	012	016


TC48L

Flamme lang • Flame long
Flamme longue

- 12 Schneiden • Blades • Lames
- 20 Schneiden • Blades • Lames
- 30 Schneiden • Blades • Lames

300.000

5



Shank	L mm		
Size	ISO	Order No.	8,0
Ø ₁₀ mm			010
FG	• 500 314 249 072...	TC48L.314...	010
	• 500 314 249 042...	TC48LF.314...	012
	• 500 314 249 032...	TC48LUF.314...	012