



## CeraPro

Stück blau-türkis, mittelgrob, HP, Figur 107, 13 mm, ISO 050



Artikelnr.: 10268  
Hersteller: Edenta  
Gewicht: 0.005 kg  
Höhe: 8 mm  
Breite: 70 mm  
Länge: 70 mm

**€ 12,82**

*Zzgl. MwSt.*

**Stück** Schleifer blau-türkis, mittelgrob, HP, Figur 107, 13 mm, ISO 050

Abrasive Trimmer with diamond particles. For fast and gentle material removal during the processing of ZrO<sub>2</sub> ceramic-/full ceramic and porcelain. For pre-contouring and finishing of larger vestibular and oral surfaces without heating the object.

Work with light pressure.

- Drehzahl 10.000 - 15.000 U/min.
- Optimale Drehzahl 12.000 U/min.
- Drehzahl (8003.150HP & 8004.120HP) 5.000 U/min.
- Verblend- und Keramiktechnik, Kronen-/ Brückentechnik, Occlusalbearbeitung Keramik